## OBSAH/CONTENT

### CZ

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## Literatura

Boeufgras, J. M., Balyer, J. M., Allard, F., Diaz, M.: A new computer program for routine interpretation of API system. *2nd conference on taxonomy and automatic identification of bacteria*. (1987) Prague, June, 29 – July, 3.

Lapage, S. P., Bascomb, S., Willcox, W. R., Curtis, M. A.: Identification of bacteria by computer. General aspects and perspectives. (1973) *J. Gen. Microbiol.* 77, 273-290

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# EN

## Introduction

Register (code book) for the kit ENTEROtest 24 N innovated is designed for the identification of bacteria by MIKRO-LA-TEST kit ENTEROtest 24 N innovated. The identification of organisms and the interpretation can be performed in three ways:

1. Manually – comparing the results with the identification table
2. Using the computer program TNW – the obtained results are automatically compared with information about the taxa. The information is stored in the database in the computer memory
3. The profile index identification. The profile index is the compressed extract of all the possible test results in a list form (Register).

All the remarked identification ways start with a database made up of several taxa in which are stored the percentage probabilities of positive results for each taxon.

A taxon can be:

- genus                *Morganella sp.*
- species             *Morganella morganii*
- subspecies        *Morganella morganii ssp. morganii*

### Profile coding

To facilitate the work it is suitable to transform the obtained results of biochemical tests into a numerical profile. The principle of coding is to condense the binary pieces of information (+/-) into a numerical profile (in a octal numeral system). The identification is based on tests of ENTEROtest 24 N innovated, supplemented with strip test for cytochromoxidase (OXI test).

At first, the tests are divided into the groups of three:

H		G		F		E		D		C	
URE		ARG		ORN		LYS		H2S		SCI	
SAL	1	SOR	2	MLB	3	CEL	4	LAC	5	TRE	6
DUL		ADO		ART		SUC		INO		RAF	

Now we have created groups of tests. Let's give a value equal to 1, 2 or 4 according to the position of the test in it's group for each test:

H		G		F		E		D		C	
URE	1	ARG	1	ORN	1	LYS	1	H2S	1	SCI	1
SAL	2	SOR	2	MLB	2	CEL	2	LAC	2	TRE	2
DUL	4	ADO	4	ART	4	SUC	4	INO	4	RAF	4

The identification is performed by summing up all the positive test values in the group and creating a number. The number has as many figures as is the count of the groups (including the incomplete group):

H		G		F		E		D		C	
URE	1 - =0	ARG	1+ =1	ORN	1 - =0	LYS	1- =0	H2S	1+=1	SCI	1+ =1
SAL	2 - =0	SOR	2- =0	MLB	2+ =2	CEL	2+=2	LAC	2+=2	TRE	2+ =2
DUL	4 + =4	ADO	4- =0	ART	4+ =4	SUC	4- =0	INO	4=0	RAF	4+ =4
	=4		=1		=6		=2		=3		=7

So the identification code is **416237**.

## Using the analytical profile index

Look for the numerical profile in the index; the profiles are arranged according to their increasing numerical values in the Register.

The profile is listed along with the following information:

- value of the two following indexes for each taxon:
  - a) **Identification score (% id.)** – an estimate saying how closely the profile corresponds to the taxon relative to all the other taxa in the database; the taxons are arranged according to their % id. values.
  - b) **T-index (Tin)** – an estimate saying how closely the profile corresponds to the most typical set of reactions for given taxon. Its value varies between 0 and 1 and is inversely proportional to the number of atypical tests.
- **The list of the tests against** the identification (T.against) for the first listed taxon – if any, followed by the percentage of positive reactions
- **The list of the supplementary tests** – if the taxa are not well distinguished by ENTEROtest 24 N innovated; the supplementary tests are listed with the percentage of positive reaction.
- **Comments** on the quality of identification are created on the base of the values of % id. and T-index:

### Identification score

% id  $\geq$  99 ... the strain is excellently distinguished

% id  $\geq$  93 ... the strain is very good distinguished

% id  $\geq$  85 ... the strain is distinguished

% id < 85 ... the strain is not distinguished

### T-index

T-index  $\geq$  0,75 ... the typical strain

T-index  $\geq$  0,50 ... the less typical strain

T-index  $\geq$  0,25 ... the atypical strain

T-index < 0,25 ... the entirely atypical strain

In the case that the **profile is not found in the analytical profile** index there may be the following reasons for this:

- The profile is too atypical and it's frequency of occurrence is very low
- The profile corresponds to the taxon, which is not included in the database

## Refernces

Boeufgras, j. M., Balyer, J. M., Allard, F., Diaz, M.: A new computer program for routine interepretation of API system. *2nd conference on taxonomy and automatic identification of bacteria*. (1987) Prague, June, 29 – July, 3.

Lapage, S. P., Bascomb, S., Willcox, W. R., Curtis, M. A.: Identification of bacteria by computer. General aspects and perspectives. (1973) *J. Gen. Microbiol.* 77, 273-290

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